



DART™

Intranasal mucosal atomisation device




Airway Management ▪ Atomisation Devices

DART™ intranasal mucosal atomisation device

High-quality mucosal atomisation in the palm of your hands

DART is indicated for atomising solutions across the naso and oropharyngeal mucous membranes and provides accurate delivery with an average particle size of approximately 50 µm. No IV line is required and solutions can be delivered quickly and painlessly. Atomised medication delivery results in therapeutic onset delivery similar to intravenous delivery. There is no special training required for a system wide adoption. The delivery method is painless for the patient, helping them remain calm in a high stress situation.

Our DART range of mucosal atomisation devices includes a number of product options to ensure you always have the optimal configuration for your requirements. This includes options both with and without a syringe and an option with a vial adaptor (a blunt plastic cannula to draw-up the medication/solution to the syringe) for convenience and ease of use. All DART devices are single patient use and packaged in boxes of 25 units.



Code	Description	Box Qty.
DART100	DART intranasal atomisation device with 3 ml syringe	25
DART110	DART intranasal atomisation device with 1 ml syringe	25
DART140	DART intranasal atomisation device with 3 ml syringe and vial adapter	25
DART300	DART intranasal atomisation device without syringe	25

Make an enquiry

DART is manufactured by Pulmodyne Inc, and distributed by Intersurgical Ltd

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS6.22 • Issue 3 10.24



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com