




O2-MAX™

Disposable CPAP device



Critical Care ▪ Pulmodyne®



Quality, innovation and choice

Interact with us




www.intersurgical.com

O2-MAX™

Continuous positive airway pressure is becoming the standard of care for patients suffering from respiratory distress. The O2-MAX System can deliver up to 140 L/min of flow to the patient, with 30% FiO₂. Customise your FiO₂ to increase from 30, 60, or 90% with the Trio adaptor. Adjustable PEEP

options includes a 3-SET™ (5/7.5/10 cmH₂O) or 5-SET™ (5/7.5/10/12.5/15 cmH₂O) to suit a variety of patient needs. Additionally, this completely disposable system offers integrated nebulisation while maintaining a consistent flow to the patient.



O2-MAX accessories

Code	Description	Box Qty.
313-8000	O2-MAX Trio control unit	10
313-8300	O2-MAX exhalation filter	10
313-7716	O2-CPAP™ fixed value PEEP valve, 5.0 cmH ₂ O, blue	10
313-7717	O2-CPAP fixed value PEEP valve, 7.5 cmH ₂ O, yellow	10
313-7718	O2-CPAP fixed value PEEP valve, 10.0 cmH ₂ O, green	10
313-7719	O2-CPAP fixed value PEEP valve, 12.5 cmH ₂ O, pink	10
313-7720	O2-CPAP fixed value PEEP valve, 15.0 cmH ₂ O, red	10
313-7721	O2-CPAP fixed value PEEP valve, 20.0 cmH ₂ O, gold	10
313-7725	3-SET adjustable PEEP valve (5/7.5/10 cmH ₂ O), white	10
313-7730	5-SET adjustable PEEP valve (5/7.5/10/12.5/15 cmH ₂ O), grey	10



O2-MAX™ Fixed System with expandable tube (delivers 30% FiO₂)

Code	Description	Box Qty.
The systems below are provided with 3-SET™ adjustable PEEP valve		
313-7555XN	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve, nebuliser and BiTrac™ ED mask, large adult	10
313-7556XN	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve, nebuliser and BiTrac ED mask, medium adult	10
313-7555X	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve and BiTrac ED mask, large adult	10
313-7556X	O2-MAX Fixed System with expandable tube, 3-SET adjustable PEEP valve and BiTrac ED mask, medium adult	10
The systems below are provided with 5-SET™ adjustable PEEP valve		
313-7554XN	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED mask, large adult	10
313-7557XN	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED mask, medium adult	10
313-7554X	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve and BiTrac ED mask, large adult	10
313-7557X	O2-MAX Fixed System with expandable tube, 5-SET adjustable PEEP valve and BiTrac ED mask, medium adult	10

O2-MAX Trio System (delivers adjustable 30%, 60%, 90% FiO₂)

Code	Description	Box Qty.
The systems below are provided with expandable tube and nebuliser		
313-8018XN	O2-MAX Trio System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED mask, large adult	10
313-8008XN	O2-MAX Trio System with expandable tube, 5-SET adjustable PEEP valve, nebuliser and BiTrac ED mask, medium adult	10
The systems below are provided with smoothbore tube		
313-8018	O2-MAX Trio System with smoothbore tube, 5-SET adjustable PEEP valve and BiTrac ED mask, large adult	10
313-8008	O2-MAX Trio System with smoothbore tube, 5-SET adjustable PEEP valve and BiTrac ED mask, medium adult	10

BiTrac ED masks with standard and CAP headgear

Code	Description	Box Qty.
313-7031	BiTrac ED mask, extra large adult	10
313-7030	BiTrac ED mask, large adult	10
313-7029	BiTrac ED mask, medium adult	10
313-7028	BiTrac ED mask, small adult	10
313-7031C	BiTrac ED mask with CAP headgear, extra large adult	10
313-7030C	BiTrac ED mask with CAP headgear, large adult	10
313-7029C	BiTrac ED mask with CAP headgear, medium adult	10
313-7028C	BiTrac ED mask with CAP headgear, small adult	10

Quik-Connect™ disposable oxygen connections



Code	Description	Box Qty.
313-111111	Quik-Connect Ohmeda gas fitting, green	30
313-121111	Quik-Connect Chemetron gas fitting, green	30
313-131111	Quik-Connect Puritan Bennett gas fitting, green	30
313-141111	Quik-Connect Schrader gas fitting, green	30
313-351111	Quik-Connect American gas fitting, white	30
313-361111	Quik-Connect DIN gas fitting, white	30
313-371111	Quik-Connect AFNOR gas fitting, white	30
313-381111	Quik-Connect BS gas fitting, white	30
313-391111	Quik-Connect UNI gas fitting, white	30

ADAPT-2-VENT™

The ADAPT-2-VENT is an adapter elbow that allows the use of the BiTrac ED mask on standard dual limb or single limb ventilator circuits. This allows the patient to move from CPAP, using the O2-MAX™ or the GO-PAP™, to other forms of ventilation without changing the patient interface.



Code	Description	Box Qty.
313-9185	ADAPT-2-VENT 30mmF non-vented elbow with anti-asphyxiation valve	10
313-9186	ADAPT-2-VENT 30mmF standard non-vented elbow	10

3-SET, 5-SET, ADAPT-2-VENT, BiTrac, O2-MAX and Quik-Connect are manufactured by Pulmodyne Inc, USA and distributed by Intersurgical Ltd
 Ohmeda is a registered trademark of GE Healthcare Finland Oy. Puritan Bennett is a registered trademark of Covidien LP.
 Chemetron is a registered trademark of Allied Healthcare Products, Inc.

Please think before you print Save energy and paper. If you must print this information sheet please print it double sided.

IS1.64 ▪ Issue 2 07.22



Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
 T: +44 (0)118 965 6300 info@intersurgical.com www.intersurgical.com

UK ▪ Ireland ▪ France ▪ Germany ▪ Spain ▪ Portugal ▪ Italy ▪ Benelux ▪ Sweden ▪ Denmark ▪ Lithuania ▪ Russia ▪ Czech Republic
 Turkey ▪ South Africa ▪ China ▪ Japan ▪ Taiwan ▪ Philippines ▪ USA ▪ Canada ▪ Colombia ▪ Australia